



Picosecond Ultrasonics (Pump/Probe)

Laser generation of acoustic waves




Generation Mechanisms

- thermoelastic
- deformation potential
- piezoelectric



The response of a material irradiated by a pulsed or chopped laser beam. Surface acoustic waves are confined to the near surface region of the sample while bulk wave travel throughout the sample.


Pump/probe technique



The pump generates an acoustic disturbance. The path delayed probe then measures acoustic strain modulation of the optical reflectivity or out-of-plane surface motion at one instant in time (instant in time ~500 fs). Time is index by increasing the optical path of the pump.


Bandwidth considerations

Conventional



- $\xi_r \sim 1\mu\text{m}$
- $\tau \sim v_r \xi_r \sim 1\text{ ns}$

Lithographic grating



- $\xi_r \sim 100\text{nm}$
- $\tau \sim v_r \xi_r \sim 50\text{ps}$

The bandwidth for bulk waves and surface acoustic waves is related to the acoustic transit time across the skin depth and optical spot size respectively.